



Collaborative Results and Information






Overview

- The goal of this presentation is to provide information on a rapid, high volume, Officially approved PRIMARY method for fat/oil extraction, that will reduce labor costs, reduce solvent usage and increase accuracy and precision.




Features and Benefits

| | |
|---|--|
| AOCS Official Approved Procedure, Am 5-04 | The system has been proven in rigorous testing |
| 90°C Extraction temperatures & Filter Bag Technology | Greatly <u>reduces</u> extraction time. |
| Totally Automatic | Free up technicians for other work |
| Automatic Solvent recycle and Recovery | Reduces solvent usage to ~2ml per sample |
| Primary Method | Can be used to monitor/calibrate secondary methods |



What Sample Types?

- **Oilseeds** – Corn, Soybean, Soy meal, Canola, Safflower, etc.
- **Meats** – Beef, Turkey, Chicken, Sausage, further processed meats, etc.
- **Animal Feeds** – Pet Foods, Cattle Feeds, Poultry Feeds, Forages, Corn Silages, etc.
- **Other Foods** – Brownie Mixes, Baked Goods, Cereals, Snack Foods, Cookies, etc.



Comparative Study

22 Sample Types – 5 reps each
 $R^2 = 0.9996$

| Sample | Conventional | | FBT | |
|------------------|--------------|------|----------|------|
| | %Fat/Oil | SD | %Fat/Oil | SD |
| Rice Hulls | 0.3 | 0.07 | 0.2 | 0.08 |
| Soybean Meal | 1.4 | 0.01 | 1.7 | 0.05 |
| Pig Starter | 1.8 | 0.05 | 1.9 | 0.11 |
| Chick Grower | 2.2 | 0.10 | 2.2 | 0.10 |
| Cattle Feed | 2.7 | 0.10 | 2.8 | 0.08 |
| Corn | 3.0 | 0.07 | 3.5 | 0.12 |
| Chicken Breast | 3.2 | 0.07 | 3.1 | 0.05 |
| Blueberry muffin | 4.6 | 0.41 | 4.7 | 0.39 |
| Oatmeal | 5.9 | 0.08 | 5.7 | 0.21 |
| Brownie Mix | 8.8 | 0.07 | 8.5 | 0.15 |
| Turkey | 8.9 | 0.11 | 8.7 | 0.07 |
| Fish Meal | 9.9 | 0.07 | 9.8 | 0.16 |
| Ham | 10.6 | 0.03 | 10.9 | 0.11 |
| Soybeans | 21.3 | 0.08 | 21.0 | 0.44 |
| Horse feed | 22.1 | 0.18 | 22.2 | 0.05 |
| Tortilla Chips | 24.2 | 0.22 | 24.2 | 0.26 |
| Ground Beef | 28.4 | 0.16 | 28.6 | 0.23 |
| Chicken Thighs | 29.1 | 0.09 | 29.2 | 0.13 |
| Sausage | 36.4 | 0.35 | 36.7 | 0.62 |
| Safflower | 40.4 | 0.22 | 39.5 | 0.20 |
| Canola | 41.4 | 0.07 | 41.7 | 0.12 |
| Cheese Curls | 43.3 | 0.06 | 43.2 | 0.29 |

AOCS Collaborative Results

Meats

Statistical Analysis provided by AOCS

| Sample Type | turkey | ham | beef ground | chicken breast | hot dog | sausage |
|---------------------------------------|--------|------|-------------|----------------|---------|---------|
| Number of laboratories | 12 | 9 | 11 | 11 | 12 | 11 |
| Number of replicates | 24 | 18 | 22 | 22 | 24 | 22 |
| Collaborative Average, Oil/Fat % | 3.2 | 11.6 | 23.8 | 2.8 | 39.5 | 25.7 |
| Certified Labs Average ^a | 3.2 | 11.3 | 23.5 | 2.7 | 39.0 | 25.0 |
| Repeatability | | | | | | |
| s(r) = repeatability std dev | 0.21 | 0.30 | 0.24 | 0.33 | 0.35 | 0.34 |
| RSD(r) = repeatability rel. std. Dev | 6.57 | 2.59 | 1.01 | 11.89 | 0.89 | 1.33 |
| r = repeatability value | 0.58 | 0.84 | 0.67 | 0.94 | 0.98 | 0.96 |
| Reproducibility | | | | | | |
| s(R) = reproducibility std dev | 0.34 | 0.30 | 0.36 | 0.33 | 0.59 | 0.51 |
| RSD(R) = reproducibility rel std. Dev | 10.84 | 2.59 | 1.49 | 11.89 | 1.49 | 1.98 |
| R = reproducibility value | 0.96 | 0.84 | 0.99 | 0.94 | 1.65 | 1.43 |

AOCS Collaborative Results

Oilseeds

Statistical Analysis provided by AOCS

| Sample Type | soybean A | canola | soybean meal | corn A | soybean B | safflower | corn B |
|---------------------------------------|-----------|--------|--------------|--------|-----------|-----------|--------|
| Number of laboratories | 11 | 9 | 12 | 12 | 11 | 9 | 12 |
| Number of replicates | 22 | 18 | 24 | 24 | 22 | 18 | 24 |
| Collaborative Average, Oil/Fat % | 20.9 | 39.0 | 1.6 | 3.3 | 19.4 | 22.5 | 3.4 |
| Certified Labs Average ^a | 21.1 | 39.7 | 1.6 | 3.6 | 19.7 | 23.0 | 3.7 |
| Repeatability | | | | | | | |
| s(r) = repeatability std dev | 0.35 | 0.23 | 0.14 | 0.31 | 0.38 | 0.53 | 0.39 |
| RSD(r) = repeatability rel. std. Dev | 1.7 | 0.6 | 8.5 | 9.5 | 1.97 | 2.36 | 11.48 |
| r = repeatability value | 0.98 | 0.65 | 0.39 | 0.88 | 1.07 | 1.49 | 1.10 |
| Reproducibility | | | | | | | |
| s(R) = reproducibility std dev | 0.63 | 0.68 | 0.27 | 0.42 | 0.62 | 0.83 | 0.41 |
| RSD(R) = reproducibility rel std. Dev | 3.0 | 1.7 | 16.3 | 12.7 | 3.19 | 3.69 | 11.93 |
| R = reproducibility value | 1.76 | 1.90 | 0.75 | 1.18 | 1.73 | 2.33 | 1.14 |

AOCS Collaborative Results

Feeds & Forages (*includes Pet Foods*)

Statistical Analysis provided by AOCS


| Sample Type | poultry starter | cattle feed | pig starter | alfalfa | cat food | dog food | corn silage |
|---------------------------------------|-----------------|-------------|-------------|---------|----------|----------|-------------|
| Number of laboratories | 11 | 10 | 11 | 11 | 12 | 12 | 12 |
| Number of replicates | 22 | 20 | 22 | 22 | 24 | 24 | 24 |
| Collaborative Average, Oil/Fat % | 3.3 | 3.2 | 5.6 | 2.4 | 6.3 | 6.8 | 2.3 |
| Certified Labs Average ^a | 3.5 | 3.0 | 5.5 | 2.2 | 6.2 | 6.9 | 2.3 |
| Repeatability | | | | | | | |
| s(r) = repeatability std dev | 0.24 | 0.18 | 0.20 | 0.39 | 0.27 | 0.35 | 0.23 |
| RSD(r) = repeatability rel. std. Dev | 7.3 | 5.6 | 3.6 | 16.1 | 4.2 | 5.23 | 9.87 |
| r = repeatability value | 0.68 | 0.51 | 0.56 | 1.08 | 0.75 | 0.99 | 0.63 |
| Reproducibility | | | | | | | |
| s(R) = reproducibility std dev | 0.42 | 0.20 | 0.28 | 0.50 | 0.30 | 0.35 | 0.51 |
| RSD(R) = reproducibility rel std. Dev | 12.6 | 6.1 | 5.0 | 20.7 | 4.7 | 5.23 | 22.45 |
| R = reproducibility value | 1.16 | 0.55 | 0.78 | 1.39 | 0.83 | 0.99 | 1.44 |

AOCS Collaborative Results


Other Foods

Statistical Analysis provided by AOCS

| Sample Type | oat meal | brownie mix | cookies | bkft cereal | tortilla chips | crackers | potato chips | cheese curls |
|---------------------------------------|----------|-------------|---------|-------------|----------------|----------|--------------|--------------|
| Number of laboratories | 12 | 12 | 11 | 12 | 12 | 10 | 11 | 12 |
| Number of replicates | 24 | 24 | 22 | 24 | 24 | 20 | 22 | 24 |
| Collaborative Average, Oil/Fat % | 5.8 | 8.7 | 22.7 | 2.3 | 19.9 | 23.8 | 32.0 | 30.6 |
| Certified Labs Average ^a | 5.7 | 8.7 | 23.1 | 2.3 | 20.0 | 24.0 | 32.0 | 30.5 |
| Repeatability | | | | | | | | |
| s(r) = repeatability std dev | 0.36 | 0.20 | 0.20 | 0.26 | 0.39 | 0.23 | 0.48 | 0.48 |
| RSD(r) = repeatability rel. std. Dev | 6.2 | 2.3 | 0.9 | 11.4 | 2.0 | 0.96 | 1.49 | 1.59 |
| r = repeatability value | 0.99 | 0.56 | 0.56 | 0.72 | 1.09 | 0.64 | 1.34 | 1.36 |
| Reproducibility | | | | | | | | |
| s(R) = reproducibility std dev | 0.54 | 0.31 | 0.20 | 0.36 | 0.48 | 0.23 | 0.52 | 0.69 |
| RSD(R) = reproducibility rel std. Dev | 9.4 | 3.5 | 0.9 | 15.7 | 2.4 | 0.96 | 1.61 | 2.27 |
| R = reproducibility value | 1.52 | 0.86 | 0.56 | 1.00 | 1.35 | 0.64 | 1.45 | 1.94 |





AOCS Collaborative *Regression Analysis*



◆ = *Certified Labs*
■ = *Collaborative Labs*

Some Key Users From Around the World



Tyson Foods, Inc.

