

Which Instrument is Right for You?

Back

XT20
Fat Analyzer

XT15
Extractor

XT10
Extractor



Instrument	XT10 Extractor	XT15 Extractor	XT20 Fat Analyzer
Extractions/run	10	15	20
Extractions/day	100	150	>200
Solvent Recovery	~ 90%	>97%	>97%
Solvent/sample	3-5ml	0.7 - 2ml	0.7 - 2ml
Automation	Semi-automatic ¹	Automatic ²	Fully Automatic
Vent System	Minimal vapor output Vent line to Hood	Charcoal Activated System No Vent Required	Charcoal Activated System No Vent Required
Controls	Microprocessor Control	Microprocessor Control	Microprocessor Control Organic Vapor Detection Fiber Optic Sensor

¹ Requires manual solvent insertion

² Requires manual lid close