


1	2	Spring	AISI 301/1.4310	600079	
2	2	Insert/back ferrule/front ferrule	AISI 316L/1.4404	600073	
3	2	Hex nut M16x1.5/M15x1.5	AISI 316L/1.4404	600077	
4	2	O-ring 7.1x1.6	Silicone	600825	
5	2	Quick coupling for teflon	AISI 316L/1.4404	800171	
6	2	O-ring 10.3x2.4	Silicone	800830	
7	1	Quick coupling with W9 hose piece	AISI 316L/1.4404	400086	
8	1 meter	Teflon tube 10/8	Teflon	900047	
POS.	Qty.	Item.	Material.	Ident no.	Remarks

keofitt
WORLD LEADERS IN STERILE SAMPLING

Keofitt a/s
Hans Egedes Vej 19
DK-5210 Odense NV
Denmark
Phone +45 63 16 70 80
Fax +45 63 16 70 81
E-mail: keofitt@keofitt.dk

Scale:	<input checked="" type="checkbox"/>	Date:	Sign.:	
1:1	<input type="checkbox"/>	Drawn	04.03.25	Is
		Sign.		
A4	<input type="checkbox"/>	Passed		

Date: Modifications Rev. Sign.



Tube for steam boiler
Complete 1 meter

Dr. No. Revision

550003