



**Where to use :** With all steam sterilizable sampling valves.  
**How to use:** Manually-Instructions available  
**Use for:** Sterilizing with steam.

**Accessories:**

**Ignition:** Piezoelectric (included)  
**Connection hose:** PTFE (included)  
**Burner:** Labogaz (included)

**Material:**

**Steel Parts:** AISI 316L(1.4404)  
**O-rings:** EPDM  
**Foot:** EPDM

**Key measurements**

<b>Drawing part A:</b>	394 mm.	15.51 In.
<b>Drawing part B:</b>	194 mm.	7.64 In.
<b>Capacity:</b>	0-250 ml.	0 – 8.5 oz.
<b>Net weight:</b>	1.8 kg./pc	3.97 lb
<b>Temp. range:</b>	1 <sup>0</sup> C - 100 <sup>0</sup> C	34°F - 212°F

Pressure and temperature can be increased by using the a Pressure relief valve (850059).

**Please note!** It is prohibited to air mail or air freight the gas cans, it is a standard product available throughout the world.

**Caution!** Open flame