

POS.	Qty.	Item.	Material.	Ident. no.	Remarks
1	1	Vessel for water, complete	AISI 316L/1.4404	550111	
10	1	O-ring 36.1x3.53	EPDM	550187	
2	1	Lever gauge, complete	-	550080	
3	1	Stand for steam boiler, complete	AISI 316L/1.4404	550001	
4	1	Ignition, complete	-	550046	
5	1	Burner stand	-	550142	t. gasconister type C206
6	1	Foot for steam boiler	EPDM	550190	
7	1	Tube for steam boiler, complete 0.5 meter	-	550002	
8	0.75 m	Chain	-	900185	
9	1	Handle for steam boiler	-	550126	

keofitt
 KEOFITTECHNIKA
 WOOD LADENS INDUSTRIELSTRASSE 10

keofitt o/a
 Koenigsgasse Weg 19
 DK-5210 Odense NV
 Phone +45 63 18 70 80
 Fax +45 63 18 70 81
 E-mail: keofitt@keofitt.dk

Steam boiler

Dwg. No.

5500000

Scale:

1:2

A3

Drawn

04.03.31

Sign.

ls

Passed


Date:

Sign.:

Remarks

Revision

Date:	Modifications	Rev.	Sign.



1 2 3 4 5 6 7 8 9 10 11 12

A B C D E F G H

1 2 3 4 5 6 7 8 9 10 11 12

A B C D E F G H