



1	2	Fitting for water gauge	AISI 316L/1.4404	550165	
2	2	O-ring 10x2	EPDM	550186	
3	2	Screw for water gauge M5x20	AISI 316L/1.4404	550169	
4	4	O-ring	EPDM	550184	
5	1	Ball for level ø6	Nylon	550163	
6	1	Glass water gauge	Polysulfon (PSU)	550159	
POS.	Qty.	Item.	Material.	Ident no.	Remarks

keofitt
 WORLD LEADERS IN STERILE SAMPLING

Keofitt a/s
 Hans Egedes Vej 19
 DK-5210 Odense NV
 Denmark
 Phone +45 63 16 70 80
 Fax +45 63 16 70 81
 E-mail: keofitt@keofitt.dk

Scale:	1:1	Date:	04.03.26	Sign.:	ls
1:1		Sign.			
A4		Passed			

Date: Modifications Rev. Sign.

Level gauge complete

Dr. No. Revision

550060

