



Data sheet:550142
Item: Labogaz C206
Burner stand for steam boiler.

Gas Cannister
 Item no.: 550145



Steam boiler
 Item no.: 550000



Where to use :	With Keofitt steam boiler	
How to use:	Insert in steam boiler.	
Material:		
Steel Parts:	Steel	
Burner:	Labogaz	
Gas canister:	C206	
Key measurements		3.54 In.
Drawing part A:	90 mm.	7.56 In.
Drawing part B:	192 mm.	0.49 lb
Net weight:	0.22 kg./pc	0.87 lb